



### Product Information

- eShield Cell Phone Size
- Reorder number: EC2100
- Overall size: 9" x 14"
- Functional size: 9" x 4.5"
- Sterilized via gamma radiation
- Contents of each peel pouch sterile unless seal is broken or package is punctured. Examine before use.
- ANSI/AAMI PB70 Level 4
- Not made with natural rubber latex
- Single use, disposable
- Do not resterilize
- Avoid direct exposure to sunlight
- Made in the USA
- 16 CFR Part 1610 Class 1 Flammability

Do you have a new product idea? Submit it to Flash!™ off our website and see your product come to life! Idea for eShield was brought to our attention by Dr. Alex Meininger. With the help of Whitney Medical Solutions Dr. Meininger was able to see his idea come to fruition.

ec2100h – Rev. 3



Whitney Medical Solutions  
847-470-9300  
www.whitneymedicalsolutions.com  
Made in U.S.A.

For questions, comments, or concerns, please contact

### Please watch video instructions prior to use




<http://qrs.ly/bd3rppq8>

For an electronic copy of this user manual in various languages, scan the following QR code



<http://qrs.ly/al3cbcz>



- Note: The entire cover is packaged sterile; "Sterile"/"Unsterile" indicator line is to help scrub person visualize which areas the circulator may have touched. Users should not have cross contact.
- Note: To maintain optical clarity of the cover, store in a cool dry place out of direct light.
- Do not use cover on devices with sharp edges that can potentially damage the cover and compromise the sterile barrier.
- If you notice a tear or other breach in the cover at any time during use, immediately discard the cover.

### General Cautions

Whitney Medical Solutions equipment drapes are to be used to cover a variety of surgical and non-surgical equipment in various settings throughout the clinical setting. These drapes are used to protect the equipment from contamination during various procedures.

### Indications for Use

Please watch video instructions prior to use




<http://qrs.ly/bd3rpq8>

## WRITTEN INSTRUCTIONS

### Cover Preparation Phase


- 1) **Scrub person:** Place hands inside cuffs and separate to create opening to cover.



 **Tip:** When putting cell phone into the cover during next step, the scrub person can hold cover horizontally to allow for circulator to easily place the cell phone inside the opening.

Note: The entire cover is packaged sterile; "Sterile"/"Unsterile" indicator line is to help scrub person visualize which areas the circulator may have touched. Users should not have cross contact.

- 2) **Circulator:** Gently place cell phone inside the cover.




Caution: Avoid dropping the cell phone into the cover. Place covered cell phone on a sterile stand for steps 3-7 so that the weight is supported.

- 3) **Scrub person:** Hold the Green *Hold Here* Circles and pull in opposite directions to tear off top. Discard top portion.



Caution: The newly created top edge of the cover is considered to be unsterile after the cuffs are torn off. Avoid contact with this area for steps 3-7.

- 4) **Scrub person:** Remove the adhesive tape liner.




- 5) **Scrub person:** Fold cover over at fold lines 1 then 2 to seal along adhesive tape. To avoid excess air, seal half of the cover, press down to remove air, and then seal other half.




Caution: Avoid unsterile areas when folding over top of cover. To help avoid cross contact, corners of cover can be folded in prior to folding over.

- 6) **Scrub person:** If "NOT SEALED" is visible, cover is not ready for use.



- 7) **Scrub person:** Once "SEALED" is formed, cover is ready for use.



- 8) **Scrub person:** Once the cover is done being used, unfold top seal to remove the cell phone. Dispose of the potentially biohazardous cover in accordance with facility protocol.