


Indications for Use

Whitney Medical Solutions equipment drapes are to be used to cover a variety of surgical and non-surgical equipment in various settings throughout the clinical setting. These drapes are used to protect the equipment from contamination during various procedures.

General Cautions

- If you notice a tear or other breach in the cover at any time during use, immediately discard the cover.
- Do not use cover on devices with sharp edges that can potentially damage the cover and compromise the sterile barrier.
- Note: To maintain optical clarity of the cover, store in a cool dry place out of direct light.
- Note: The entire cover is packaged sterile; "Sterile"/"Unsterile" indicator line is to help scrub person visualize which areas the circulator may have touched. Users should not have cross contact.



<http://qrs.ly/al3cbc>



For an electronic copy of this user manual in various languages, scan the following QR code



Please watch video instructions prior to use

<http://qrs.ly/1q3rpbqz>

Whitney Medical Solutions
6153C W. Mulford St.
Niles, IL 60714
847-470-9300
www.whitneymedicalsolutions.com



Do you have a new product idea? Submit it to Flash!™ off our website and see your product come to life! Idea for eShield was brought to our attention by Dr. Alex Meininger. With the help of Whitney Medical Solutions Dr. Meininger was able to see his idea come to fruition.

ec2400h – Rev. 1



Product Information

- eShield SLR Camera Size
- Reorder number: EC2400
- Overall size: 14" x 21"
- Functional size: 14" x 10.5"
- Sterilized via gamma radiation
- Contents of each peel pouch sterile unless seal is broken or package is punctured. Examine before use.
- ANSI/AAMI PB70 Level 4
- Not made with natural rubber latex
- Single use, disposable
- Do not resterilize
- Avoid direct exposure to sunlight
- Made in the USA
- 16 CFR Part 1610 Class 1 Flammability
- Sterile rubber band attached to peel pouch

Please watch video instructions prior to use




<http://qrs.ly/1q3rpqz>

WRITTEN INSTRUCTIONS:

Camera Preparation Phase


1) Circulator: Prepare camera for use; remove camera strap, lens cap, extend lens barrel to max extension, if applicable, extend flash and orient display screen.



Cover Preparation Phase


2) Scrub person: Place hands inside cuffs and separate to create opening to cover.



 Tip: When putting camera into the cover during next step, the scrub person can hold cover horizontally to allow for circulator to easily place the camera inside the opening.

Note: The entire cover is packaged sterile; "Sterile"/"Unsterile" indicator line is to help scrub person visualize which areas the circulator may have touched. Users should not have cross contact.

3) Circulator: Gently place camera inside the cover. Place camera so that camera is sitting flat on the bottom and lens is facing the indicated lens side.




Caution: Avoid dropping the camera into the cover. Place covered camera on a sterile stand for steps 4-11 so that the weight is supported.

4) Scrub person: Hold the Green Hold Here Circles and pull in opposite directions to tear off top. Discard top portion.




Caution: The newly created top edge of the cover is considered to be unsterile after the cuffs are torn off. Avoid contact with this area for steps 5-8.

5) Scrub person: Remove the adhesive tape liner




6) Scrub person: Fold cover over at fold lines 1 then 2 to seal along adhesive tape. To avoid excess air, seal half of the cover, press down to remove air, and then seal other half.




Caution: Avoid unsterile areas when folding over top of cover. To help avoid cross contact, corners of cover can be folded in prior to folding over.

7) Scrub person: If "NOT SEALED" is visible, cover is not correctly sealed.




8) Scrub person: Once "SEALED" is formed, cover is correctly sealed.





9) Circulator: Open Sterile Rubber Band peel pouch and remove rubber band.




 Tip: An extra rubber band is included in case one needs to be discarded.


10) Scrub person: Lay camera flat with lens facing up and hold cover tight against lens to place rubber band around barrel.



 Tip: Have SLR lens extended as far out as possible to allow for easier manipulation once the rubber band is applied.

11) Scrub person: Cover is now ready for use.



 Tip: Extra rubber band can be used to secure any excess material.

12) Scrub person: Once the cover is done being used, carefully remove rubber band from camera bezel and unfold top seal to remove the camera. Dispose of the potentially biohazardous cover in accordance with facility protocol.