

Iso-Tray™

Draw Blood in Isolation Rooms Safely



Phlebotomy In Isolation Rooms

Iso-Tray™ is a disposable phlebotomy tray for carrying blood drawing tools and accessories into isolation rooms where contamination risks are high.

Features & Benefits

- 10 dedicated spots for blood sample tubes.
- 2 dedicated spots for blood culture tubes.
- Iso-Trays store flat - requires less space and has interlocking flaps for easy assembly.
- Iso-Trays are made from sturdy paperboard, a renewable resource.
- Center well for all other blood drawing supplies.
- Aseptic Non-touch Technique compatible.
- Reduces the risk of infection during invasive clinical procedures.



Disposable



Single-use procedure tray promotes asepsis per ANTT*

Economical



Practical



Carries everything you need and eliminates extra trips!

*Aseptic Non-Touch Technique