



Manage Your Biohazards with Whitney Medical Solutions:  
Real Solutions for Real Problems™



## Collect and Dispose of Lab Waste Easily, and Safely:

Sturdy corrugated paperboard construction offers strength, and the 2 mil red plastic liner is leak-resistant.

Dispose of non-sharp biohazardous lab waste, including petri dishes, tubes, plates, pipets and pipet tips, specimen cups, and tubing sets.



### CAUTION:

Not a sharps container—do not use for disposal of contaminated sharps.

### Biohazard Keeper™ Floor Model:

Floor Model	
Order No.:	BH2009
Dimensions:	12"W x 12"D x 27"H
Capacity:	64 Liters
Shipping Weight:	14 Lbs.
Quantity:	Case of 6
<b>Description:</b> Sturdy corrugated construction, heavy-duty 2 mil liner, and unique wide-mouth opening with Reclose-A-Lid / Remove-A-Lid™ cover.	

### Lid Features:

- Unique Reclose-A-Lid / Remove-A-Lid™ reduces odor and waste visibility
- Removable, wide-mouth opening for large items
- A separate top flap permanently closes the container for final disposal



### Additional Features:

- Heavy-duty 2 mil liner
- Sturdy corrugated construction
- Stores flat
- Assembles in seconds
- Clearly marked biohazard symbol
- Warning labels in English and Spanish



You'll wonder how you ever managed *without us.*™