



Which eShield™ Cover is Right for My Application?

		
Purpose:	To keep handheld electronics clean and reduce cross-contamination.	To allow non-sterile handheld electronics to be used safely in a sterile field.
Applications and Environments for Use:	Clean / Disinfected Environments (partial list): hospitals, dental offices, schools, medical offices, sales and corporate meeting interactions, hotels, rental equipment	Sterile Environments (partial list): hospitals, operating rooms, industrial clean rooms, research facilities, tissue / sample processing, pharmaceutical manufacturing and testing
Devices:	Smartphones, tablets, and small electronics (remote controls, small keyboards, etc.)	Smartphones, tablets, digital cameras, SLR cameras, and small electronics (handheld medical or laboratory devices)
Sterile:	No Packaged non-sterile, 100 covers per carton.	Yes Individually packaged in sterile packages, 25 per carton.
Certifications and Standards:	FDA Approved Materials	FDA Approved Class II Medical Device Bio-compatibility Test: Non-Cytotoxic ASTM F1670 Certified Meets FDA and Industry Standards for Leak Resistance
Assembly:	Single-User Assembly 2-person technique is not required.	2-Person Assembly Unique tear-off cuff design allows proper 2-person aseptic technique during set-up. Double-fold closure says "Sealed" when correctly folded
Clarity:	Clear Film Provides clear view of screen.	Ultra-Clear Film Ultra clear polyethylene film for visibility and clarity for video and photography applications.
Touch Screen Sensitivity:	Yes Screen is responsive to touch with finger and while wearing single or double latex gloves. Repositionable tape allows adjustment for a tight fit and better touch screen performance.	Yes Screen is responsive to touch while wearing single or double latex gloves.
Additional Photography Features:	No	Yes Ultra-clear material allows for excellent photography. Double-faced adhesive rings hold lens bezel to wall of the cover for minimal internal reflection.
Learn More:	Visit www.StayCleanCovers.com Call 800-338-4237, or email us at info@whitneym.com	Visit www.SterileCover.com Call 800-338-4237, or email us at info@whitneym.com



Manufactured by Whitney Medical Solutions,
The Experts in Germ-Free Digital.

