

BIOHAZARD PIPET KEEPER™ CONTAINERS



Dispose of pipets, swabs, and other long biohazardous objects easily

PIPET KEEPER™

More economical, holds items up to 15" long



PIPET KEEPER™ XL

75% more capacity and holds items up to 16" long



Size Comparison

Specifications:

MADE IN THE USA

Order No.:	BH2006
Dimensions:	3"W x 3"H x 15½"L
Capacity:	2300 cc
Quantity:	Case of 50
Shipping Weight:	8 Lbs.

Specifications:

MADE IN THE USA

Order No.:	BH2016
Dimensions:	3.5"W x 4.5"H x 16½"L
Capacity:	4200 cc
Quantity:	Case of 50
Shipping Weight:	13 Lbs.

SHARED FEATURES & BENEFITS

- **Sturdy Safety:** Strong paperboard, a renewable resource, printed with orange biohazard symbols
- **Plastic Lined:** Avoid soak-through with webbed corners to hold moderate amounts of liquid
- **Disposable:** Save time and make it easier to maintain high lab standards
- **Priced Lower:** Use instead of rigid plastic sharps containers
- **Interlocking, Reclosable Top:** Easy to use whether you are right or left-handed
- **Shipping & Storage Efficiency:** Ships flat and can easily be folded together for use



Real Solutions for Real Problems